

BEFORE / AFTER | 2 months post 1 tx | courtesy of Yonick Plastic Surgery



# BBL<sup>®</sup> + Halo<sup>®</sup> + SkinTyte<sup>™</sup>

BBL: 560 nm filter, 10 J/cm<sup>2</sup>, 10 ms, 15 °C Cooling, 2 Passes; 560 nm filter, 17 J/cm<sup>2</sup>, 15 ms, 20 °C Cooling, 2 Passes

515 nm filter, 12 J/cm<sup>2</sup>, 10 ms, 20 °C Cooling, 2 Passes, 11mm round adaptor

SkinTyte: 590 nm filter, 10 W/cm<sup>2</sup>, 12 seconds, 30 °C Cooling

Halo: 1470 nm: 325 μm, 40%; 2940 nm: 70 μm, 32%

